

HDS Microscope Slides for Sysmex SP-1000i



Product Features

- Designed specifically for use with Sysmex SP-1000i
- Manufactured in Germany from premium quality extra white glass
- Precleaned, non-coated 100% wettable surface for consistent quality of blood smears
- Thermal printer paint optimised for use with ribbon printer in SP-1000i
- Clipped corners to prevent jamming during slide transportation
- Long side ground edge of 90 ° and a short side ground edge of 45 ° to avoid jamming into the integrated CellaVision imaging magazine
- Certified to ISO 13485 (medical devices quality management system)
- CE mark certified for IVD (In-Vitro Diagnostics) medical devices
- FDA registered

Specifications

- Slide size - 25.4 mm x 75.5 mm
- Slide thickness - 1.0 mm
- Long side ground edge 90 °
- Short side ground edge 45 °
- Clipped corners 2.0mm
- 50 slides per box
- 1,500 per carton (30 boxes)

Code	Description
Slides for Automated Preparation	
HD 61296SM 1P0	HDS Slides, Colourfrost, White, for Sysmex SP-1000i

Contact your local Grale HDS office for technical support, product advice or to request product samples for evaluation.

► Sydney **1800 257 213** | Melbourne **03 9874 8577** | enquiries@grale.com.au | www.grale.com.au