

## A Pathology Supply Partner Like No Other

### Fresh New look for Grale HDS Slides and Coverslips



Grale HDS is proud to introduce a new packaging design for our premium range of microscope slides and coverslips. The outer packaging may have changed, however the quality of the slides and coverslips have not. The source of supply has not changed and you will continue to receive a premium quality product manufactured in Germany.

- Clear brand identification – providing quick recognition of the trusted Grale HDS brand.
- Product type – clearly identified on top of the box flap.
- Product batch details – detailing batch/lot number and expiry date
- Product code – for easy reordering when you are running low.
- Coverslips are supplied in a custom foam within a plastic box with a hinged lid and firm clip front closure.

Product specifications, box quantities and prices are not changing, so you can continue to order the same quality products using the current product codes.

The fresh new look Grale HDS slide and coverslip packaging will be phased in from late 2014.



Contact your local Grale HDS office for technical support, product advice or to request product samples for evaluation.

► Sydney **1800 257 213** | Melbourne **03 9874 8577**


## Grale HDS Microscope Slides



- Manufactured in Germany from premium quality extra white glass.
- Double frosted or Colourfrost, with plain, 45 ° or 90 ° ground edges.
- Colourfrost slides are available with a marking area in six colours.
- Slides are cellophane wrapped as standard and packed in vacuum sealed cartons.
- Slide size - 76 mm x 26 mm.
- 50 slides per box.
- 1,500 slides per carton (30 boxes).

Code	Description
<b>Plain Edge Slides</b>	
HD 21102 1P0	Grale HDS Slides, Plain Edge, Double Frosted
HD 21206 1P0	Grale HDS Slides, Colourfrost, Plain Edge, White
HD 21207 1P0	Grale HDS Slides, Colourfrost, Plain Edge, Blue
HD 21208 1P0	Grale HDS Slides, Colourfrost, Plain Edge, Pink
HD 21209 1P0	Grale HDS Slides, Colourfrost, Plain Edge, Yellow
HD 21210 1P0	Grale HDS Slides, Colourfrost, Plain Edge, Green
HD 21211 1P0	Grale HDS Slides, Colourfrost, Plain Edge, Orange
<b>45 ° Ground Edge Slides</b>	
HD 41104 1P0	Grale HDS Slides, Ground Edge 45 °, Double Frosted
HD 41206 1P0	Grale HDS Slides, Colourfrost, Ground Edge 45 °, White
HD 41207 1P0	Grale HDS Slides, Colourfrost, Ground Edge 45 °, Blue
HD 41208 1P0	Grale HDS Slides, Colourfrost, Ground Edge 45 °, Pink
HD 41209 1P0	Grale HDS Slides, Colourfrost, Ground Edge 45 °, Yellow
HD 41210 1P0	Grale HDS Slides, Colourfrost, Ground Edge 45 °, Green
HD 41211 1P0	Grale HDS Slides, Colourfrost, Ground Edge 45 °, Orange
<b>90 ° Ground Edge Slides</b>	
HD 41808 1P0	Grale HDS Slides, Ground Edge 90 °, Double Frosted
HD 41906 1P0	Grale HDS Slides, Colourfrost, Ground Edge 90 °, White
HD 41907 1P0	Grale HDS Slides, Colourfrost, Ground Edge 90 °, Blue
HD 41908 1P0	Grale HDS Slides, Colourfrost, Ground Edge 90 °, Pink
HD 41909 1P0	Grale HDS Slides, Colourfrost, Ground Edge 90 °, Yellow
HD 41910 1P0	Grale HDS Slides, Colourfrost, Ground Edge 90 °, Green
HD 41911 1P0	Grale HDS Slides, Colourfrost, Ground Edge 90 °, Orange
<b>Adhesive Slides</b>	
HD 41300 1P0	Grale HDS Slides, Colourfrost, Ground Edge 90 °, Adhesive
HD 41400 1P0	Grale HDS Slides, Colourfrost, Ground Edge 90 °, Polysine
HD 41800 1P0	Grale HDS Slides, Colourfrost, Ground Edge 90 °, Advanced Adhesive
<b>Special Application Slides</b>	
HD 61296SM 1P0	Grale HDS Slides, Colourfrost, Clipped, White
HD BC9045CC 1P0	Grale HDS Slides, Ground Edge 90 °, Bevelled End, Clipped, Thermal Customized
HD 41906 1P0C	Grale HDS Slides, Colourfrost, Ground Edge 90 °, Clipped, White

## Grale HDS Microscope Coverslips



- Manufactured in Germany, available in standard No. 1 thickness, in vacuum sealed cartons.
- Suitable for use in automatic coverslippers.

Code	Description
HD LD2222 1.01P0	Grale HDS Coverslips, No1, 22 x 22 mm, 2000 in Tray, 10 x Pk200
HD LD2232 1.01P0	Grale HDS Coverslips, No1, 22 x 32 mm, 1000 in Tray, 10 x Pk100
HD LD2240 1.01P0	Grale HDS Coverslips, No1, 22 x 40 mm, 1000 in Tray, 10 x Pk100
HD LD2250 1.01P0	Grale HDS Coverslips, No1, 22 x 50 mm, 1000 in Tray, 10 x Pk100
HD LD2440 1.01P0	Grale HDS Coverslips, No1, 24 x 40 mm, 1000 in Tray, 10 x Pk100
HD LD2450 1.01P0	Grale HDS Coverslips, No1, 24 x 50 mm, 1000 in Tray, 10 x Pk100
HD LD2455 1.01P0	Grale HDS Coverslips, No1, 24 x 55 mm, 1000 in Tray, 10 x Pk100
HD LD2460 1.01P0	Grale HDS Coverslips, No1, 24 x 60 mm, 1000 in Tray, 10 x Pk100